**12** RESPONSIBLE  
CONSUMPTION  
AND PRODUCTION

Vestfrost Solutions is working towards reaching the UN - Global Sustainable Development Goals by 2030.

The Sustainable Development Goals are the blueprint to achieve a better and more sustainable future for all.

In order to implement Goal no 12 "Responsible Consumption and Production", this manual has been printed on recycled paper.

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# Warning

As the appliance contains flammable refrigerant, as stated on nameplate, it is essential to ensure that the refrigerant pipes are not damaged.

**The quantity and type of the refrigerant used in your appliance is indicated on the rating plate.**

Standard EN378 specifies that the room in which you install your appliance must have a volume of 1m<sup>3</sup> per 8 g of hydro-carbon refrigerant used in the appliances. This is to avoid the formation of flammable gas/air mixtures in the room where the appliance is located in the event of a leak in the refrigerant circuit.

## **WARNING:**

Ventilation openings in the appliance or in built-in structures must be kept clear.

## **WARNING:**

Do not use other mechanical devices or means to accelerate the defrosting process or to remove rime other than those recommended by the manufacturer.

## **WARNING:**

Do not damage the refrigerant system.

## **WARNING:**

The appliance may not be exposed to rain.

## **WARNING:**

This appliance is not intended for use by persons (including children) with reduced physical, sensory or mental capabilities, or lack experience and knowledge, unless they have been given supervision or instructions concerning use of the appliance by a person responsible for their safety. Children should be supervised to ensure that they do not play with the appliance

## **WARNING:**

Children must not play with, on, or around the appliance.

## **WARNING:**

Children must not clean the appliance or carry out general maintenance unless they are at least 8 years old and are being supervised.

**WARNING:**

Always, keep the keys in a separate place and out of reach of children.

**WARNING:**

Before servicing or cleaning the appliance, switch off circuit breaker.

**WARNING:**

Cleaning and user maintenance shall not be made by children without supervision

**WARNING:**

Danger risk of fire or explosion.  
Flammable refrigerant used,  
as stated on nameplate. To  
be repaired only by trained  
personnel.

**WARNING:**

Do not store explosive substances such as aerosol cans with a flammable propellant in this appliance.

**WARNING:**

When positioning the appliance, ensure the supply cord is not trapped or damaged.

**WARNING:**

Do not locate multiple portable socket-outlets or portable power supplies at the rear of the appliance.

**Contents**

Warning .....	2
Introduction.....	4
Operational guidelines.....	4
Get to know your VFS 048 SDD .....	5
Freezing waterpacks .....	5
Placing and unpacking .....	6
Installation and start up .....	7
Control panel .....	8
Maintenance and cleaning.....	9
Defrosting .....	10
Trouble-shooting.....	11
Warranty, spare parts and service .....	12
Wiring diagram .....	13
Spare parts .....	14
Disposal.....	16

## Introduction

The philosophy behind the VFS 048 SDD system is:

- Simplicity
- Safety
- Robustness
- Reliability

The VFS 048 SDD system is composed of:

- VFS 048 SDD freezer
- Solar panels
- Installation equipment for the solar panels
- Wires from the solar panels to the VFS 048 SDD

The VFS 048 SDD is a waterpack freezer. It has one compartment for freezing and storage of waterpacks.

The freezing capacity is 3x 0,6 liter waterpacks pr 24h at 43°C ambient. The storage capacity is 29 waterpacks.

With correct installation and by following a few maintenance routines, the system will work trouble-free for many years.

## Operational guidelines

The compressor requires that the solar panels are able to deliver a power above a certain threshold in order to start. When started, the compressor is able to operate at lower power consumption than the starting threshold. The compressor may need a couple of starting attempts to get running. So be patient!

During days with changing weather conditions, the compressor may start and stop several times.

The noise coming from the VFS 048 SDD appliance is low, and primarily comes from a fan cooling the compressor.

The compressor is almost noiseless. The fan starts a moment before the compressor and stops a moment after the compressor.

The controller is factory set and should be adjusted only by authorized personnel.

# Get to know your VFS 048 SDD


fig. 1

## Freezing waterpacks

Place waterpacks in the slots – after 24h the packs are frozen

Frozen waterpacks can be stored in the middle of the freezer


Fig. 2

# Placing and unpacking

## Unpacking

Unpack the appliance and check that it has not been damaged. If you observe any damage, please inform your supervisor.

## Placing

Install the appliance in a dry and well ventilated place. Avoid installation near heat sources or in direct sunlight.

Place the appliance on the floor and make sure that it is level.

### Important:

**Do not use open fire or cigarettes in the proximity of the appliance.**

## Adjustable feet

Level the appliance by screwing the adjustable feet up or down.

If the appliance is to be placed on a soft surface, e.g. floorboards or a carpet, it is best to recheck whether the appliance is still level after a period of time as the underlying surface may give under the weight of the appliance.


fig. 3

**The Climate Class is indicated on the rating plate.**

### Climate class:

0, 1, 2, 3, 4, 6 ,8 = ambient 32°C

### Climate class

5 and 7 = ambient 43°C


fig. 4

## Ventilation

It is important that the appliance is placed well ventilated, so that air can circulate unhindered above, below, and around the appliance. Fig. 2 illustrates how the necessary air circulation can be ensured. There must be at least 30 mm clearance between the base of the appliance and the floor.

Ventilation openings in the appliance or in built-in structures must be kept clear.

**NOTE: There must be free access to the ventilation grille. (On the backside.)**

# Installation and start up

Connect the appliance to the solar panels according to table 1, which states the lengths and sizes of the wires used for connecting the appliances.

Standard supplied cable size and length will be 6mm<sup>2</sup> and 20m.

The screened cables coming from the PV panels are fitted with a male/female MC4 compatible connectors. The connections are then to be mounted to the prefixed male/female MC4 compatible connectors placed at the back of the appliance - See fig. 4

The connections are of plug and play type and can only be fitted correct - male to female/female to male.

Clean the appliance thoroughly before start up.

Size:		Max. length:	
Cross section [mm <sup>2</sup> ]	AWG [Gauge]	[m]	[ft]
2,5	12	8	26
4	12	13	42
6	10	20	65
10	8	32	105

Table 1

## Start up

Before the appliance is loaded with water-packs, the ice pack ballast must be frozen (see page 5 fig. 2). The cool-down time in an ambient temperature of 43°C will be approx. 5-7 days.


Fig. 5

# Control panel


fig. 6

The green diode indicates that there is power.

# Solar panels

## Placing

Solar panels are always to be placed on a stable roof or bar, which is safe and easily accessible (not for the public).

The solar panels are always to be fitted according to the supplier's instructions and with the screws and fittings included or recommended by the manufacturer.

Do not place the solar panels in the shadows of buildings, trees, etc. The surfaces of the solar panels are to be turned directly towards the sun.

## Specifications

The specifications which the solar panels are to be observed are available in table 2.

<b>Type :</b>	Type 2 according E003/PV01 direkt drive
<b>Nominal voltage:</b>	18V
<b>Nominal power:</b>	360W (standard test conditions)
<b>Number of parallel connected modules per array:</b>	4
<b>Warranty:</b>	Min. 2 years
<b>Performance guarantee:</b>	25 years (to 80% of original value)
<b>Certificates:</b>	UL, TÜV or compliance with the requirements of IEC 61215

Table 2

# Maintenance and cleaning

## Daily maintenance:

Check lid is properly closed and fits tightly to the cabinet of the appliance.

## Monthly maintenance:

Clean the grille on the back side of the refrigerator once each month.

## Yearly maintenance:

Electric connections and components are to be checked and cleaned once a year or more if necessary.

## Maintenance of the solar panels

### Weekly maintenance:

The solar panels are to be cleaned for dust once each week – or as required. Cleaning must be carried out with water using a soft cloth or wash rag.

### Monthly maintenance:

Avoid shading

It should be regularly monitored that shades do not come up such as new towering trees, as this will decrease the amount of energy produced by the system.

### Yearly maintenance:

Electric connections and components are to be checked and cleaned at least once a year or more often if required. Because of risk of fire, it is necessary to remove dust and dirt. Checking that PV panels are clean, free of fractures, scratches, corrosion, moisture penetration and browning.

Cabling should be checked to ensure it is secure.

Checking the mounting hardware to ensure it is in good condition and ensuring the earth connection is continuous.

Checking of junction boxes to ensure there is no water accumulation and that the integrity of lid seals, connections and clamping devices is intact.

All bolts and nuts are to be re-tightened 2 weeks after installation, and then once a year.

## Cleaning

Disconnect the power supply before cleaning by using the on/off switch placed in the front of the appliance. The best way to clean the appliance is by using luke warm water with a small amount of unscented detergent. Never use cleaning agents that scour. Use a soft cloth. Rinse with clean water and dry thoroughly. It is important to prevent water from running into the control panel.

The sealing strip around the lid must be cleaned regularly to prevent discolouration and prolong service life. Use clean water. After cleaning the sealing strip, check that it continues to provide a tight seal.

If the appliance is not being used for any period of time, switch the appliance off, disconnect the power supply, empty the appliance, clean the inside, and leave the lid open to allow air circulation and prevent smells.

# Defrosting

The water-packs stored and the humidity of the air gradually causes frost and ice to form in the freezer compartment. A thin layer of frost or ice does not affect freezer performance, but a thick layer could cause the freezer not to cool water-packs satisfactorily. Small amounts of loose frost can be removed using a plastic scraper.

Never use electrical apparatus or sharp implements as these may damage the inside of the appliance.

When the layer of frost and ice exceeds 4-5 mm the freezer should be completely defrosted.

Defrosting should be performed in the morning between 7am to 9am, just before the solar radiation again is powerful enough to start compressor and run refrigeration system.

Proceed as follows:

1. Disconnect the power supply by switching off the ON/OFF button. Important! ONLY switch off the power supply of the freezer compartment. (See fig. 7)
2. Remove all water-packs from freezer compartment. Place them in an empty coldbox or wrap them up to keep them as cold as possible.
3. Keep the freezer lid open.
4. Open the two drainage plugs. One is place inside the freezer compartment, and one placed in left side front of appliance.
5. Place a tray under the defrost water drain.
6. Place a bowl of hot (but not boiling) water in the freezer, close the lid and wait approx. 30 minutes before removing loose frost with the plastic scraper.

7. To prevent water from damaging the floor during defrosting, it is advisable to check the defrost water tray from time to time and to cover the floor around the defrost water drain to protect against ice and splashes. When all the frost and ice have melted, clean the inside cabinet, as well lid, lid gasket, and outside of the cabinet using water with a mild, perfume-free detergent. Dry the freezer with a soft cloth, and make sure no moisture is remaining after cleaning.
8. Reinsert the internal and external drainage plugs, and check lid gasket is not faulty.
9. Reconnect the power supply by the ON/OFF button.
10. Check waterpacks are not leaking before returning them to the freezer compartment. Place waterpacks according instructions, approx. 5 minutes after having switched it on. Close the lid.


fig. 7

# Trouble-shooting

Fault	Possible cause	Remedy
Compressor is not running, and the ice packs are not cold	Be patient, it is most likely that the compressor will start within a few minutes.	If this is not the case, check the following: - Check that power is connected and that the wire from the solar panel to the appliance is intact. - Check the fuse and replace it if necessary. - If the above is OK, call technical supervisor.
Compressor is running, and the temperature is too high	The ventilation grill is blocked.  The lid is not closed properly.  The fan in the compressor compartment is blocked or defective.  The temperature in the room in which the appliance is installed is too high.	Ensure unhindered air circulation.  Ensure that the lid is closed properly.  Check that the fan is running, if not it should be replaced. See section: Service  Shield the appliance against direct sun light and ensure more ventilation to the room.
No temperature is displayed	There is not enough light for the solar sensor.	Turn on the light.
No light in green diode during day time	Switch is turned OFF  Fuse is burned  Diode is defect	Turn on the switch  Replace the fuse.  Replace diode.

Table 3

# Warranty, spare parts and service

## Warranty disclaimer

Faults and damage caused directly or indirectly by incorrect operation, misuse, insufficient maintenance, incorrect building, installation or mains connection. Fire, accident, lightening, voltage variation or other electrical interference, including defective fuses or faults in mains installations are not covered by the warranty.

Repairs performed by others than approved service centres and any other faults and damage that the manufacturer can substantiate are caused by reasons other than manufacturing or material faults are not covered by the warranty.

Please note that changes to the construction of the appliance or changes to the component equipment of the appliance will invalidate warranty and product liability, and the appliance cannot be used lawfully. The approval stated on rating plate will also be invalidated.

Transport damage discovered by the buyer is primarily a matter to be settled between the buyer and the distributor, i.e. the distributor must ensure that such complaints are resolved to the buyer's satisfaction.

Before calling for technical assistance, please check whether you are able to rectify the fault yourself. If your request for assistance is unwarranted, e.g. if the appliance has failed as a result of a blown fuse or incorrect operation, you will be charged the costs incurred by your call for technical assistance.

## Spare parts

When ordering spare parts, please state the type, serial and product numbers of your appliance. This information is given on the rating plate. The rating plate contains various technical information, including type and serial numbers.




Fig. 8

# Wiring diagram

# VFS 048 SDD

Position	Item number	Item name
0061	0-6538001	Filter drier,
0087	0-6038175	Base plate fittings, complete
0129	8-036510214	Compressor - BD35K 101Z0211 10-45V
2170	0-A9301260103	Distance piece
5000	3010049	Drain plug
5013	3040400	Adjustable foot
5052	2042052-01	Mounting plate for condenser VLS024
5057	6010436	WOT condenser
5144	5000917024	Lid foamed without handle/hinges
5402	1510059	Hinge
5450	3011135-01	Top part for hinge cover
5451	3010032-01	Bottom part for hinge cover
5570	2042044	Motor screen
5605	304090501	Handle with lock (push and turn)
5606	8090342-94	Inlay for handle,
5682	1510046	Key, set - 2pc.
5683	2040145	Catch for handle
5684	3010265-01	Cover for catch
5685	600098801	Lock with keys,(push and turn)
5709	7060104	Frame for cover rocker switch bezel
5710	6520038	Thermostat, freezer
5718	7060105	Rocker Switch Cover
5727	6520845	Electronics for BD35K compressor
5730	7020245	"on-off" switch
5739	7010139	Status indicator
5755	7060062	Fuse holder
5756	7020382-03	Thermometer, solar
5770	3010005-01	Thermostat knob
5801	8470156	Fuse 15A
5863	A903007	"Well-nut", Screw,
6221	7030823	Fan

# Disposal

## *Information for Users on Collection and Disposal of Old Equipment and used Batteries*

 These symbols on the products, packaging, and/or accompanying documents mean that used electrical and electronic products and batteries should not be mixed with general household waste. For proper treatment, recovery and recycling of old products and used batteries, please take them to applicable collection points, in accordance with your national legislation and the Directives 2012/19/EU and 2006/66/EC.

By disposing of these products and batteries correctly, you will help to save valuable resources and prevent any potential negative effects on human health and the environment which could otherwise arise from inappropriate waste handling.

For more information about collection and recycling of old products and batteries, please contact your local municipality, your waste disposal service or the point of sale where you purchased the items.

Penalties may be applicable for incorrect disposal of this waste, in accordance with national legislation.

### **For business users in the European Union.**

If you wish to discard electrical and electronic equipment, please contact your dealer or supplier for further information.

### **Information on Disposal in other Countries outside the European Union**

These symbols are only valid in the European Union. If you wish to discard this product, please contact your local authorities or dealer and ask for the correct method of disposal.

### **Note for the battery symbol:**

 This symbol might be used in combination with a chemical symbol. In this case it complies with the requirement set by the Directive for the chemical involved.

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