


VESTFROST
SOLUTIONS

SERVICEMANUAL - EN

MODEL - ULTF-C74i ULTF-C198i ULTF-C296i ULTF-C383i



BIOMEDICAL SOLUTIONS

FREEZERS

CONTENTS

ATTENTION!	4
MAINTENANCE SCHEDULE	4
REPLACEMENT OF PARTS	5
Display - with integrated keyboard	5
Handle and inlay	6
Hinge and lid	8
Fan, Battery, Transformer, Power board	10
ULTF-C198i, ULTF-C296i, ULTF-C383i.	10
Fan, Battery, Transformer, Electronic control unit ULTF-C74i	12
COOLING SYSTEM	14
Change compressor	14
Recharging refrigerant from a service bottle.	15
ELECTRICAL SYSTEM	16
Wire diagram ULTF-C198i ULTF-C296i, ULTF-C383i	16
Wire diagram ULTF-C74i.	17
SPARE PARTS	18
Spare parts exploded view ULTF-C74i	18
Spare parts list ULTF-C74i	19
Spare parts exploded view ULTF-C198i	20
Spare parts list ULTF-C198i	21
Spare parts exploded view ULTF-C296i	22
Spare parts list ULTF-C296i	23
Spare parts exploded view ULTF-C383i	24
Spare parts list ULTF-C383i	25
LIST OF PARAMETERS	26
List of Ultra low parameters	26
List of LCP parameters	27
List of LCS parameters	28

ATTENTION!

Before servicing or cleaning the appliance, unplug the appliance from the mains or disconnect the electrical power supply.

ATTENTION!

This service manual is intended for use by persons having electrical and mechanical training and a level of knowledge of these subjects generally considered acceptable in the appliance repair trade. The manufacturer cannot be responsible,

nor assume any liability, for injury or damage of any kind arising from the use of this manual.

ATTENTION!

Bare pcbs are sensitive to ESD, so please take precautions to ensure that tools, persons and appliance is equalised.

ATTENTION!

Compressor and tubes can be hot, please allow the unit to cool down before repairing.


MAINTENANCE SCHEDULE

All below maintenance is recommended.
Please investigate the local legislation in your area.


TASK	SCHEDULE
Cleaning of interior	Weekly
Cleaning of compressor compartment (including condenser and compressor)	Monthly
Cleaning of sealing in the lid	Quarterly
Replacement of battery for back up	Every third year

REPLACEMENT OF PARTS

Display - with integrated keyboard




Release the display by pushing in the plast lock (located under display) and lifting up the panel simultanously




Carefully demount the plug/wire from the display back side

Mount display in reverse order


Handle and inlay




Use a small flat-blade screwdriver



Tilt out the inlay




Use the Torx 20 screwdriver and remove the screw on the right side of the handle




Turn the lever up (1)

Pull the aluminum tube to the right (2)




The right plastic part and the aluminum tube
are now removed




Use the Torx 20 screwdriver and remove the
screw on the left side of the handle

Mount the new handle in reverse order


Hinge and lid



Remove hinge covers




Pull out the bottom hinge cover




Pull up on the top hinge cover


Important The lid must be fully open before proceeding




Use the Torx 20 screwdriver to loosen 4 screws in each hinge.



The lid is removed with a light lift.



Use the Torx 20 screwdriver to remove 4 screws in each hinge




Mount the new hinde/lid in reverse order

IMPORTANT check that the lid flush with the cabinet (mounted straight) ...

... that the sealing strip closes tightly. Check-with a piece of paper


Fan, Battery, Transformer, Power board

ULTF-C198i, ULLT-C296i, ULLT-C383i




For access to the engine room compartment, remove engine room cover. Engine room cover is fixed with seven built-in clips

2 pcs at each side and 3 at the bottom




To release the clip function use a small flat-blade screwdriver

For each released clips, gently pull out the cover.




The battery is removed by removing the cable tie with eg cutter plier, and disconnect the red and black wire.


W




Remove the nuts and wire connectors and the box with transformer and power board holder can be pulled out of the compressor compartment



Remove the power board by removing the 5 screws (one in the middle and one in each corner of the power board). The transformer is removed from the holder by removing 2 screws.




Unplug fan from power board and Remove the 2 bolts with hexagon head




To remove the fan from the bracket, remove 2 bolts with hexagon head


Fan, Battery, Transformer, Electronic control unit ULTF-C74i




For access to the engine room compartment, remove engine room cover. Engine room cover is fixed with 2 screws




Unplug fan from power board and Remove the 2 bolts with hexagon head



The battery is removed by removing the bracket (2 screws), and disconnect the red and black wire.




Remove 2 nuts and wire connectors and the box with transformer and power board holder can be pulled out of the compressor compartment



Remove the power board by removing the 5 screws (one in the middle and one in each corner of the pore board). The transformer is removed from the holder by removing 2 screws.

COOLING SYSTEM

Change compressor



Repairs to the cooling system may only be carried out by trained professionals.

When changing the compressor, the filter drier must also be replaced.


It is important to cut out the filter drier, as heating with solder flame will release moisture inside the dry filter. This will vaporize and drive moisture from filter drier back into system.

The quantity and type of the refrigerant used in your appliance is indicated on the rating plate.

Do not make sharp bending angles on the pipes.

The distance between the pipes must be at least 12 millimetres.

Recharging refrigerant from a service bottle




The recharging shall take place in a room where the ambient temperature is between 20 and 25 °C.

1. Connect service bottle, compressor and vacuum pump (A) to a T
2. Open valve (C) and let vacuum pump run for minimum 30 min
3. Close valve (C)
4. Place the bottle on a scale if possible.
5. Open valve (B) slowly
6. When the refrigerant stops flowing into the system, start the compressor.
7. Let the compressor run for 5 min
8. Check amount of refrigerant charged on the scale, should be approx. the same as on the rating plate
9. Repeat from 5 if necessary


Note! If you start the compressor at point 8 with the valve open then you can ruin the compressor.

ELECTRICAL SYSTEM

Wire diagram ULLT-C198i ULLT-C296i, ULLT-C383i




Wire diagram ULTF-C74i



SPARE PARTS


Spare parts exploded view ULTF-C74i



Spare parts list ULTF-C74i

Position	Item number	Item name
0061	0-6538001	Filter drier,
0087	0-6038175	Base plate fittings, complete
0129	8-036510297	Compressor - Complete NEU2155U Embraco
0898	7090048	Sensor
3011	8471472	Display
3104	7090466	Transformer
5054	2040229	Fittings for condenser
5055	6010135	Condenser
5144	5200908640002	Lid foamed without handle/hinges
5145	6000986	Lid, lower380 x 380 x 27
5346	7095666	Electronic control unit
5399	8471798	Wheels with break
5402	1510059	Hinge
5450	3011135-01	Top part for hinge cover
5451	3010032-01	Bottom part for hinge cover
5548	7-049954504	Nature-R - refilling with filter
5562	6520467	Terminal Board Complet
5605	304090501	Handle with lock (push and turn)
5606	8090314	Inlay for handle,
5682	1510231	Key, set - 2pc.
5683	2040145	Catch for handle
5684	3010265-01	Cover for catch
5685	600098801	Lock with keys,(push and turn)
5860	7060046	Battery
5875	6520867	Dry contact for Alarm complete / Terminal box 3 pol
5882	6520393	Cover
5972	2042987	Motor screen
6553	7535084	Mains lead
7100	7090478	Fan


Spare parts exploded view ULTF-C198i



Spare parts list ULTF-C198i

Position	Item number	Item name
0061	0-6538001	Filter drier,
0087	0-6038175	Base plate fittings, complete
0126	8-036510213	Compressor -SC15CNX.2
0898	7090048	Sensor
3011	8471472	Display
3104	7090466	Transformer
5000	3010049	Drain plug
5054	2040467	Fittings for condenser
5055	6010135	Condenser
5112	5000905206651	Lid foamed without handle/hinges
5346	7095667	Electronic control unit
5399	8471797	Wheels with break
5402	1510135	Hinge with spring
5404	5020211	Foam lid/Inner lid 463x300x29,5 mm
5450	3011135-01	Top part for hinge cover
5451	3010032-01	Bottom part for hinge cover
5533	A921052	Electrical parts for compressor SC15CNX.2
5548	7-049954505	Nature-R - refilling with filter
5605	304090501	Handle with lock (push and turn)
5606	8090054	Inlay for handle,
5682	1510238	Nikel-plated Key
5683	2040145	Catch for handle
5684	3010265-01	Cover for catch
5685	600098801	Lock with keys.(push and turn)
5740	5020287	Lid for ice bank
5860	7060046	Battery
5875	6520867	Dry contact for Alarm complete / Terminal box 3 pol
5970	3010308-01	Motor screen
5971	3520001	Motor screen
6553	7535084	Mains lead
7100	7090478	Fan


Spare parts exploded view ULTF-C296i



Spare parts list ULTF-C296i

Position	Item number	Item name
0061	0-6538001	Filter drier,
0087	0-6038175	Base plate fittings, complete
0129	8-036510002	Compressor complete SC21CNX.2 - 104H8166
0898	7090048	Sensor
3011	8471472	Display
3104	7090466	Transformer
5000	3010049	Drain plug
5052	2040504	Fittings for condenser
5053	6010068	WOT condenser
5118	500090465101	Lid foamed without handle/hinges
5346	7095668	Electronic control unit
5399	8471797	Wheels with break
5402	1510135	Hinge with spring
5405	5020212	Foam lid/Inner lid 463x363x29,5 mm
5450	3011135-01	Top part for hinge cover
5451	3010032-01	Bottom part for hinge cover
5533	A921053	Electrical parts for compressor SC15CNX.2
5548	7-049954508	Nature-R - refilling with filter
5605	304090501	Handle with lock (push and turn)
5606	8090054	Inlay for handle,
5682	1510238	Nikel-plated Key
5683	2040145	Catch for handle
5684	3010265-01	Cover for catch
5685	600098801	Lock with keys,(push and turn)
5860	7060046	Battery
5875	6520867	Dry contact for Alarm complete / Terminal box 3 pol
5970	3010308-01	Motor screen
5971	3520001	Motor screen
6553	7535084	Mains lead
7100	7090478	Fan

Spare parts exploded view ULTF-C383i



Spare parts list ULTF-C383i

Position	Item number	Item name
0061	0-6538001	Filter drier,
0087	0-6038175	Base plate fittings, complete
0129	8-036510002	Compressor complete SC21CNX.2 - 104H8166
0898	7090048	Sensor
3011	8471472	Display
3104	7090466	Transformer
5000	3010049	Drain plug
5052	2040504	Fittings for condenser
5053	6010068	WOT condenser
5122	500090365101	Lid foamed without handle/hinges
5346	7095669	Electronic control unit
5399	8471797	Wheels with break
5402	1510136	Hinge with spring
5404	5020211	Foam lid/Inner lid 463x300x29,5 mm
5405	5020212	Foam lid/Inner lid 463x363x29,5 mm
5450	3011135-01	Top part for hinge cover
5451	3010032-01	Bottom part for hinge cover
5533	A921053	Electrical parts for compressor SC15CNX.2
5548	7-049954509	Nature-R - refilling with filter
5605	304090501	Handle with lock (push and turn)
5606	8090054	Inlay for handle,
5682	1510238	Nikel-plated Key
5683	2040145	Catch for handle
5684	3010265-01	Cover for catch
5685	600098801	Lock with keys,(push and turn)
5860	7060046	Battery
5875	6520867	Dry contact for Alarm complete / Terminal box 3 pol
5970	3010308-01	Motor screen
5971	3520001	Motor screen
6553	7535084	Mains lead
7100	7090478	Fan

LIST OF PARAMETERS

List of Ultra low parameters

Profile name	PROFILE 1	Ultra low			
Description	Temperature range: -86C to -20C Set point: -82C				
INDEX	PAR NAME	VALUE	UNIT		
1	ADR	***	***	49	AD1
2	ALL	50	Tenth of°C/F	50	AD2
3	ALH	50	Tenth of°C/F	51	AD3
4	AAL	-999	Tenth of°C/F	52	AD4
5	AAH	-500	Tenth of°C/F	53	CDT
6	ALD	***	***	54	FCN
7	ADD	***	***	55	FCX
8	DTX	***	***	56	EDT
9	DTD	***	***	57	FEN
10	BOD	***	***	58	FEX
11	VOP	***	***	59	HPO
12	ISP	***	***	60	HUM
13	OSP	***	***	61	HRH
14	DIN	***	***	62	HRL
15	DOO	***	***	63	VBH
16	OF1	***	***	64	VBL
17	OF2	***	***	65	VBX
18	SET	-820	Tenth of°C/F	66	VBY
19	HYH	***	***	67	CSO
20	HYL	***	***	68	CSH
21	SLL	-860	Tenth of°C/F	69	CSD
22	SLH	-200	Tenth of°C/F	70	CNS
23	CCD	***	***	71	CXS
24	ADS	***	***	72	DCS
25	ADL	***	***	73	FSC
26	CON	***	***	74	FSA
27	COF	***	***	75	FSB
28	CPH	***	***	76	TPD
29	COP	***	***	77	TMI
30	CSP	***	***	78	HYA
31	TCO	***	***	79	OFS
32	BCO	***	***	80	LBO
33	PTD	***	***	81	LBC
34	FCE	***	***	82	LFU
35	MCT	***	***	83	LFD
36	FOP	***	***	84	LOD
37	FAS	500	Tenth of°C/F	85	PO1
38	FAD	***	***	86	PO2
39	FSD	***	***	87	PO3
40	LET	***	***	88	PI1
41	DOP	***	***	89	PI2
42	DCM	***	***	90	PO1
43	DRP	***	***	91	OP2
44	ITD	***	***	92	OP3
45	DTE	***	***	93	OP4
46	DTO	***	***	94	OP5
47	DEO	***	***		
48	P30	***	***		

List of LCP parameters

INDEX	PAR NAME	VALUE	UNIT		
1	ADR	1	Unit	48 DDT	-250 Tenth of?C/F
2	ALL	50	Tenth of?C/F	49 DDP	30 Min
3	ALH	50	Tenth of?C/F	50 AD1	24 Time hour
4	AAL	-999	Tenth of?C/F	51 AD2	24 Time hour
5	AAH	999	Tenth of?C/F	52 AD3	24 Time hour
6	ALD	10	Min	53 CDT	100 Tenth of?C/F
7	ADD	120	Min	54 FCN	30 %
8	DTX	200	Tenth of?C/F	55 FCX	100 %
9	DTD	10	Min	56 EDT	100 Tenth of?C/F
10	BOD	10	Min	57 FEN	35 %
11	VOP	0	Bit flag	58 FEX	100 %
12	ISP	129	Bit flag	59 HPO	0 %
13	OSP	33	Bit flag	60 HUM	20 %
14	DIN	0	Unit	61 HRH	0 %
15	DOO	120	Sec	62 HRL	0 %
16	OF1	0	Tenth of?C/F	63 VBH	14700 mV
17	OF2	0	Tenth of?C/F	64 VBL	12800 mV
18	SET	-820	Tenth of?C/F	65 IBC	50 mA
19	HYH	10	Tenth of?C/F	66 VBY	0 Unit
20	HYL	0	Tenth of?C/F	67 CSO	60 Min
21	SLL	-860	Tenth of?C/F	68 CSH	10 Min
22	SLH	-200	Tenth of?C/F	69 CSD	10 Min
23	CCD	5	Sec	70 CNS	66 Hz
24	ADS	15	Sec	71 CXS	150 Hz
25	ADL	1	Min	72 DCS	3 Hz
26	CON	5	Min	73 FSC	50 Tenth of?C/F
27	COF	3	Min	74 FSA	3 Thousandth of unit
28	CPH	100	%	75 FSB	20 Tenth of unit
29	COP	0	Bit flag	76 TPD	30 Sec
30	CSP	-550	Tenth of?C/F	77 TMI	-250 Tenth of?C/F
31	TCO	10000	Millisec	78 HYA	0 Tenth of?C/F
32	BCO	3	Sec	79 OFS	0 Tenth of?C/F
33	PTD	0	Tenth of?C/F	80 LBO	50 %
34	FCE	350	Tenth of?C/F	81 LBC	10 %
35	MCT	580	Tenth of?C/F	82 LFU	1000 Millisec
36	FOP	544	Bit flag	83 LFD	1000 Millisec
37	FAS	500	Tenth of?C/F	84 LOD	5 Sec
38	FAD	60	Sec	85 PO1	0 Unit
39	FSD	500	Tenth of?C/F	86 PO2	0 Unit
40	LET	-250	Tenth of?C/F	87 PO3	0 Unit
41	DOP	0	Bit flag	88 PI1	0 Unit
42	DCM	0	Min	89 PI2	0 Unit
43	DRP	60	Sec	90 OP1	0 Tenth of?C/F
44	ITD	8	Hour	91 OP2	0 Tenth of?C/F
45	DTE	500	Tenth of?C/F	92 OP3	0 Tenth of?C/F
46	DTO	30	Min	93 OP4	0 Tenth of?C/F
47	DEO	0	Min	94 SPX	0 Tenth of?C/F

List of LCS parameters

INDEX	PAR NAME	VALUE	UNIT
1	PCT	-1000	Tenth of°C/F
2	PHT	500	Tenth of°C/F
3	PCH	50	Tenth of°C/F
4	PHH	50	Tenth of°C/F
5	PCD	10	Sec
6	PHD	10	Sec
7	PCN	5	Unit
8	PCW	60	Min
9	PHN	5	Unit
10	PHW	60	Min
11	CS1	1	Bit flag
12	CS2	0	Bit flag
13	CS3	0	Bit flag
14	OR1	0	Tenth of°C/F
15	OR2	0	Tenth of°C/F
16	OP1	0	Tenth of°C/F
17	OP2	0	Tenth of°C/F
18	OP3	0	Tenth of°C/F
19	OP4	0	Tenth of°C/F
20	OP5	0	%
21	OP6	0	Tenth of°C/F
22	OP7	0	Tenth of°C/F
23	OP8	0	Tenth of°C/F
24	OP9	0	Tenth of°C/F
25	SAP	2	Unit



Revision date: 19-10-2022
Drawing no.: 8195121 rev 01