

Installation Qualification	Qualification Documentation Art. no.	VESTFROST <i>SOLUTIONS</i>
ULTF S. No.		

Installation Qualification

Test Protocol for Qualification

ULTF 500 / ULTF 700

- ☐ First qualification
- ☐ Requalification

Aim of IQ:

Practical test if the unit is installed according to the installation guidelines of the manufacturer.

Carried out by:	Date:
Checked by:	Date:

Generation date:	Date of alteration:
Author	

Installation Qualification	Qualification Documentation Art. no.	VESTFROST SOLUTIONS
ULTF S. No.		

Content

IQ 1.	DOCUMENTATION APPROVAL	3
IQ 1.1	Prior to performing the IQ	3
IQ 1.2	After performing the IQ	3
IQ 2.	PERSONNEL INVOLVED WITH COMPLETING THE IQ TESTS	4
IQ 3.	LIST OF DOCUMENTS IN APPENDIX	5
IQ 4.	UNIT IDENTIFICATION	6
IQ 4.1	Unit description	6
IQ 4.2	Markings on the unit	6
IQ 4.3	Identification acc. to type plate	7
IQ 5.	ACCESSORIES AND OPTIONAL EQUIPMENT	11
IQ 5.1	Standard accessories	11
IQ 5.2	Documentation supplied with the qualification folder	12
IQ 6.	LOCATION OF THE UNIT	13
IQ 7.	RESPONSIBLE PERSON FOR THIS UNIT	14
IQ 8.	USER'S DOCUMENTS	15
IQ 8.1	Operation log book	15
IQ 8.2	Standard operation procedures (SOP)	16
IQ 9.	INSTALLATION	17
IQ 10.	CONNECTION TO THE ZERO-VOLTAGE RELAY ALARM CONTACT	19
IQ 11.	ELECTRICAL POWER SUPPLY	21
IQ 12.	RESULT OF THE IQ-TEST PROTOCOL	22
IQ 12.1	Summary assessment and comments	22
IQ 12.2	Summary of IQ test results	23

Carried out by:	Date:
Checked by:	Date:

Generation date:	Date of alteration:
Author	

Installation Qualification	Qualification Documentation Art. no.	VESTFROST <i>SOLUTIONS</i>
ULTF S. No.		

IQ 1. Documentation approval

IQ 1.1 Prior to performing the IQ

Checked / Approved	Name / Department / Company	Date	Signature

IQ 1.2 After performing the IQ

Checked / Approved	Name / Department / Company	Date	Signature

Carried out by:	Date:
Checked by:	Date:

Generation date:	Date of alteration:
Author	

Installation Qualification	Qualification Documentation Art. no.	VESTFROST <i>SOLUTIONS</i>
ULTF S. No.		

IQ 2. Personnel involved with completing the IQ tests

The following table serves for clear identification of signatures/initials of all persons involved with the IQ tests.

Name	Initials	Department / Company	Date	Signature

Carried out by:	Date:
Checked by:	Date:

Generation date:	Date of alteration:
Author	

Installation Qualification	Qualification Documentation Art. no.	VESTFROST <i>SOLUTIONS</i>
ULTF S. No.		

IQ 3. List of documents in appendix

Description	Location	Appendix no.

Comments:

Carried out by:	Date:
Checked by:	Date:

Generation date:	Date of alteration:
Author	

Installation Qualification	Qualification Documentation Art. no.	VESTFROST <i>SOLUTIONS</i>
ULTF S. No.		

IQ 4. Unit identification

Aim of test:

The following tests serve to check the identity of the unit.

IQ 4.1 Unit description

Ultra-low temperature freezer 352/528 Eco Premium with microprocessor temperature control for long-term storage of samples in the ultra-low temperature range.

IQ 4.2 Markings on the unit

	Existant and easily legible	
	Yes	No
Type plate	<input type="checkbox"/>	<input type="checkbox"/>
CE conformity marking	<input type="checkbox"/>	<input type="checkbox"/>

Comments:

Carried out by:	Date:
Checked by:	Date:

Generation date:	Date of alteration:
Author	

Installation Qualification	Qualification Documentation Art. no.	VESTFROST SOLUTIONS
ULTF S. No.		

IQ 4.3 Identification acc. to type plate

ULTF 500

		Yes	No
Model	ULTF 500	<input type="checkbox"/>	<input type="checkbox"/>
#		<input type="checkbox"/>	<input type="checkbox"/>
Serial No.		<input type="checkbox"/>	<input type="checkbox"/>
Inventory No.			
Year of construction		<input type="checkbox"/>	<input type="checkbox"/>
Nominal temperature	-90 °C -130 °F	<input type="checkbox"/>	<input type="checkbox"/>
IP protection	IP 20	<input type="checkbox"/>	<input type="checkbox"/>
Nominal power	1.6 kW	<input type="checkbox"/>	<input type="checkbox"/>
Nominal voltage (+/-10%) at indicated power frequency	230 V / 50 Hz	<input type="checkbox"/>	<input type="checkbox"/>
Current type	1 N ~	<input type="checkbox"/>	<input type="checkbox"/>
Nominal current	7.0 A	<input type="checkbox"/>	<input type="checkbox"/>
Max operating pressure in the refrigerating system	28 bar	<input type="checkbox"/>	<input type="checkbox"/>
Filling weight of refrigerant stage 1 R290	0.15 kg	<input type="checkbox"/>	<input type="checkbox"/>
Filling weight of refrigerant stage 2 R170	0.15 kg	<input type="checkbox"/>	<input type="checkbox"/>

Comments:

Carried out by:	Date:
Checked by:	Date:

Generation date:	Date of alteration:
Author	

Installation Qualification	Qualification Documentation Art. no.	VESTFROST <i>SOLUTIONS</i>
ULTF S. No.		

ULTF 700

		Yes	No
Model	ULTF 700	<input type="checkbox"/>	<input type="checkbox"/>
#		<input type="checkbox"/>	<input type="checkbox"/>
Serial No.		<input type="checkbox"/>	<input type="checkbox"/>
Inventory No.			
Year of construction		<input type="checkbox"/>	<input type="checkbox"/>
Nominal temperature	-90 °C -130 °F	<input type="checkbox"/>	<input type="checkbox"/>
IP protection	IP 20	<input type="checkbox"/>	<input type="checkbox"/>
Nominal power	1.6 kW	<input type="checkbox"/>	<input type="checkbox"/>
Nominal voltage (+/-10%) at indicated power frequency	230 V / 50 Hz	<input type="checkbox"/>	<input type="checkbox"/>
Current type	1 N ~	<input type="checkbox"/>	<input type="checkbox"/>
Nominal current	7.0 A	<input type="checkbox"/>	<input type="checkbox"/>
Max operating pressure in the refrigerating system	28 bar	<input type="checkbox"/>	<input type="checkbox"/>
Filling weight of refrigerant stage 1 R290	0.15 kg	<input type="checkbox"/>	<input type="checkbox"/>
Filling weight of refrigerant stage 2 R170	0.15 kg	<input type="checkbox"/>	<input type="checkbox"/>

Comments:

Carried out by:	Date:
Checked by:	Date:

Generation date:	Date of alteration:
Author	

Installation Qualification	Qualification Documentation Art. no.	VESTFROST <i>SOLUTIONS</i>
ULTF S. No.		

ULTF 500 (UL Model 120V)

		Yes	No
Model	ULTF 500 (120V)	<input type="checkbox"/>	<input type="checkbox"/>
#		<input type="checkbox"/>	<input type="checkbox"/>
Serial No.		<input type="checkbox"/>	<input type="checkbox"/>
Inventory No.			
Year of construction		<input type="checkbox"/>	<input type="checkbox"/>
Nominal temperature	-90 °C -130 °F	<input type="checkbox"/>	<input type="checkbox"/>
IP protection	IP 20	<input type="checkbox"/>	<input type="checkbox"/>
Nominal power	1.4 kW	<input type="checkbox"/>	<input type="checkbox"/>
Nominal voltage (+/-10%) at indicated power frequency	115 V / 60 Hz	<input type="checkbox"/>	<input type="checkbox"/>
Current type	1 N ~	<input type="checkbox"/>	<input type="checkbox"/>
Nominal current	11.7 A	<input type="checkbox"/>	<input type="checkbox"/>
Max operating pressure in the refrigerating system	28 bar	<input type="checkbox"/>	<input type="checkbox"/>
Filling weight of refrigerant stage 1 R290	0.15 kg	<input type="checkbox"/>	<input type="checkbox"/>
Filling weight of refrigerant stage 2 R170	0.15 kg	<input type="checkbox"/>	<input type="checkbox"/>

Comments:

Carried out by:	Date:
Checked by:	Date:

Generation date:	Date of alteration:
Author	

Installation Qualification	Qualification Documentation Art. no.	VESTFROST SOLUTIONS
ULTF S. No.		

ULTF 700 (UL Model 120V)

		Yes	No
Model	ULTF 700	<input type="checkbox"/>	<input type="checkbox"/>
#		<input type="checkbox"/>	<input type="checkbox"/>
Serial No.		<input type="checkbox"/>	<input type="checkbox"/>
Inventory No.			
Year of construction		<input type="checkbox"/>	<input type="checkbox"/>
Nominal temperature	-90 °C -130 °F	<input type="checkbox"/>	<input type="checkbox"/>
IP protection	IP 20	<input type="checkbox"/>	<input type="checkbox"/>
Nominal power	1.4 kW	<input type="checkbox"/>	<input type="checkbox"/>
Nominal voltage (+/-10%) at indicated power frequency	115 V / 60 Hz	<input type="checkbox"/>	<input type="checkbox"/>
Current type	1 N ~	<input type="checkbox"/>	<input type="checkbox"/>
Nominal current	11.7 A	<input type="checkbox"/>	<input type="checkbox"/>
Max operating pressure in the refrigerating system	28 bar	<input type="checkbox"/>	<input type="checkbox"/>
Filling weight of refrigerant stage 1 R290	0.15 kg	<input type="checkbox"/>	<input type="checkbox"/>
Filling weight of refrigerant stage 2 R170	0.15 kg	<input type="checkbox"/>	<input type="checkbox"/>

Comments:

Carried out by:	Date:
Checked by:	Date:

Generation date:	Date of alteration:
Author	

Installation Qualification	Qualification Documentation Art. no.	VESTFROST <i>SOLUTIONS</i>
ULTF S. No.		

IQ 5. Accessories and optional equipment

Aim of test:

The following test serves to check the scope of delivery of the unit and the accompanying material and documentation. It is checked if the material complies with the given data and if the documentation is complete.

IQ 5.1 Standard accessories

	Supplied	Missing
Short user manual with information where to download the instruction manual	<input type="checkbox"/>	<input type="checkbox"/>
Ethernet interface	<input type="checkbox"/>	<input type="checkbox"/>
USB interface	<input type="checkbox"/>	<input type="checkbox"/>
3 shelves and 12 shelf holders with 6 screws	<input type="checkbox"/>	<input type="checkbox"/>
DIN plug for the zero-voltage relay alarm output (connected)	<input type="checkbox"/>	<input type="checkbox"/>
Set of 2 spacers for rear wall distance.	<input type="checkbox"/>	<input type="checkbox"/>

	Yes	No
Accessories completely present	<input type="checkbox"/>	<input type="checkbox"/>

If not:

Measures taken:

Carried out by:	Date:
Checked by:	Date:

Generation date:	Date of alteration:
Author	

Installation Qualification	Qualification Documentation Art. no.	VESTFROST <i>SOLUTIONS</i>
ULTF S. No.		

IQ 5.2 Documentation supplied with the qualification folder

	Supplied	
	Yes	No
Certificate DIN EN ISO 9001	<input type="checkbox"/>	<input type="checkbox"/>
EU Conformity declaration of the manufacturer	<input type="checkbox"/>	<input type="checkbox"/>
Test protocol for Installation Qualification IQ	<input type="checkbox"/>	<input type="checkbox"/>
Test protocol for Operational Qualification OQ	<input type="checkbox"/>	<input type="checkbox"/>
Wiring diagram and refrigeration cycle diagram with list of components	<input type="checkbox"/>	<input type="checkbox"/>
Spare parts list	<input type="checkbox"/>	<input type="checkbox"/>
Service manual	<input type="checkbox"/>	<input type="checkbox"/>
Instruction manual	<input type="checkbox"/>	<input type="checkbox"/>

	Yes	No
Documentation completely present	<input type="checkbox"/>	<input type="checkbox"/>

If not:

Measures taken:

Carried out by:	Date:
Checked by:	Date:

Generation date:	Date of alteration:
Author	

Installation Qualification	Qualification Documentation Art. no.	VESTFROST <i>SOLUTIONS</i>
ULTF S. No.		

IQ 6. Location of the unit

Aim of test:

This chapter defines the site of installation of the unit. The location must be the same when executing the qualification.

Note: Generally, no qualification can be assigned to other locations!

Address	
Building	
Department	
Room No.	
User group/Laboratory area	
User	

Comments:

Carried out by:	Date:
Checked by:	Date:

Generation date:	Date of alteration:
Author	

Installation Qualification	Qualification Documentation Art. no.	VESTFROST <i>SOLUTIONS</i>
ULTF S. No.		

IQ 8. User's documents

Aim of test:

This chapter lists if an operating book or one or more SOP were applied and defines their location. Further criteria (e.g., access right) can be defined. These determinations are optional.

IQ 8.1 Operation log book

	Yes	No
An operation log book was applied.	<input type="checkbox"/>	<input type="checkbox"/>

If yes:

Identification number	
Designation	
Depository	

Comments:

Carried out by:	Date:
Checked by:	Date:

Generation date:	Date of alteration:
Author	

Installation Qualification	Qualification Documentation Art. no.	VESTFROST <i>SOLUTIONS</i>
ULTF S. No.		

IQ 8.2 Standard operation procedures (SOP)

	Yes	No
One or several standard operation procedures) (SOP) were applied.	<input type="checkbox"/>	<input type="checkbox"/>

If yes:

	SOP	SOP (optional)
Identification, version umber		
Period of validity or expiration date		
Depository		

	SOP (optional)	SOP (optional)
Identification, version umber		
Period of validity or expiration date		
Depository		

Comments:

Carried out by:	Date:
Checked by:	Date:

Generation date:	Date of alteration:
Author	

Installation Qualification	Qualification Documentation Art. no.	VESTFROST SOLUTIONS
ULTF S. No.		

IQ 9. Installation

Aim of test:

The following tests serve to assure correct installation of the unit.

See instruction manual chap. 4 and 5

	ok	Not ok	Corrected
Transport seals of the shelves and any transportation protection devices and adhesives inside and outside the unit and at the doors were removed.	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
The unit is visually clean and undamaged	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
The operating manuals and accessory equipment were removed from the interior	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
The freezer is set up inside a building on an even surface, free from vibration and in a well-ventilated, dry location.	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
Blocking of the front castors with the breaks	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
Plane installation (control by means of a spirit level).	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
Spacers for rear wall distance were installed	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
Minimum distance of 100 mm / 3.94 in to the ventilation openings on the chamber front and rear observed	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
Lateral minimum distance to the next unit 250 mm / 9.84 in	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
Lateral minimum wall distance laterally, on the side without door hinge: 100 mm / 3.94 in	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
Lateral minimum wall distance laterally, on the side with door hinge: 240 mm / 9.45 in	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
Ambient humidity 70% r.H. max. non-condensing Measured value: _____ % r.H.	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
Max. ambient temperature 32 °C / 89.6 °F Measured value: _____ °C	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>

Carried out by:	Date:
Checked by:	Date:

Generation date:	Date of alteration:
Author	

Installation Qualification	Qualification Documentation Art. no.	VESTFROST SOLUTIONS
ULTF S. No.		


IQ 10. Connection to the zero-voltage relay alarm contact

Aim of test:

The following tests serve to assure correct connection to the zero-voltage relay alarm contact.

See chap. 15.5 of the instruction manual

The zero-voltage relay alarm output switches when the red LED Δ lights up on the controller display and in case of a power failure.



Description of the central alarm facility

The alarm contacts shall not be used at the moment.	<input type="checkbox"/>
--	--------------------------

In case the alarm contacts are used:

The transfer and evaluation of the alarm contacts should be described following:


Description:

Carried out by:	Date:
Checked by:	Date:

Generation date:	Date of alteration:
Author	

Installation Qualification	Qualification Documentation Art. no.	VESTFROST <i>SOLUTIONS</i>
ULTF S. No.		

Sketch of location and wiring:



Comments:

Carried out by:	Date:
Checked by:	Date:

Generation date:	Date of alteration:
Author	

Installation Qualification	Qualification Documentation Art. no.	VESTFROST <i>SOLUTIONS</i>
ULTF S. No.		

IQ 12.2 Summary of IQ test results

	Designation of test / title	Test passed		Date / Name (abbrev.)
		Yes	No	
IQ 1.	Documentation approval	<input type="checkbox"/>	<input type="checkbox"/>	
IQ 1.1.	Prior to performing the IQ	<input type="checkbox"/>	<input type="checkbox"/>	
IQ 1.2.	After performing the IQ	<input type="checkbox"/>	<input type="checkbox"/>	
IQ 2.	Personnel involved with completing the IQ tests	<input type="checkbox"/>	<input type="checkbox"/>	
IQ 3.	List of documents in appendix	<input type="checkbox"/>	<input type="checkbox"/>	
IQ 4.	Unit identification	<input type="checkbox"/>	<input type="checkbox"/>	
IQ 4.1	Unit description	<input type="checkbox"/>	<input type="checkbox"/>	
IQ 4.2	Markings on the unit	<input type="checkbox"/>	<input type="checkbox"/>	
IQ 4.3	Identification acc. to type plate	<input type="checkbox"/>	<input type="checkbox"/>	
IQ 5.	Accessories and optional equipment	<input type="checkbox"/>	<input type="checkbox"/>	
IQ 5.1	Standard accessories	<input type="checkbox"/>	<input type="checkbox"/>	
IQ 5.2	Documentation supplied with the qualification folder	<input type="checkbox"/>	<input type="checkbox"/>	
IQ 6.	Location of the unit	<input type="checkbox"/>	<input type="checkbox"/>	
IQ 7.	Responsible person for this unit	<input type="checkbox"/>	<input type="checkbox"/>	

Carried out by:	Date:
Checked by:	Date:

Generation date:	Date of alteration:
Author	

Installation Qualification	Qualification Documentation Art. no.	VESTFROST <i>SOLUTIONS</i>
ULTF S. No.		

	Designation of test / title	Test passed		Date / Name (abbrev.)
		Yes	No	
IQ 8.	User's documents	<input type="checkbox"/>	<input type="checkbox"/>	
IQ 8.1	Operation log book	<input type="checkbox"/>	<input type="checkbox"/>	
IQ 8.2	Standard operation procedures (SOP)	<input type="checkbox"/>	<input type="checkbox"/>	
IQ 9.	Installation	<input type="checkbox"/>	<input type="checkbox"/>	
IQ 10.	Connection to the zero-voltage relay alarm contact	<input type="checkbox"/>	<input type="checkbox"/>	
IQ 11.	Electrical power supply	<input type="checkbox"/>	<input type="checkbox"/>	
IQ 12.	Result of the IQ-Test protocol	<input type="checkbox"/>	<input type="checkbox"/>	
IQ 12.1	Summary assessment and comments	<input type="checkbox"/>	<input type="checkbox"/>	
IQ 12.2	Summary of IQ test results	<input type="checkbox"/>	<input type="checkbox"/>	

Carried out by:	Date:
Checked by:	Date:

Generation date:	Date of alteration:
Author	