

+20°C

BETWEEN
-86°C

BIOMEDICAL DEVICES

**COMPACT
SERIES**

PERFECT IN-LAB SOLUTION



Small footprint



Compact series from +20°C to -86°C

Many laboratories have limited space and the new Compact series feature personal freezers and refrigerators which fit perfectly on the counter at the work station. In developing the Compact series, it has been a point of emphasis to ensure that the models have low noise and are user friendly, ensuring a health work environment as well as being easy to use in the daily work.

True to our commitment towards a sustainable future, the Compact series has been made to have no negative impact on the environment and are marked with our "Environmental friendly – Icon" as your guarantee.

Larger sized cabinets on the market uses more energy, but if the capacity they offer is not really necessary it is worth considering exploring the smaller Compact series from Vestfrost Solutions. Energy savings can be upwards of EUR 60 annually.

COMPACT SERIES



Specifications



R-90

SPECIFICATIONS

	R-90
Dimensions, H x W x D, mm	652 x 595 x 642
Gross volume, Litres	90
Net volume, Litres	76
Temperature range, C°	+2 to +20
Number of shelves	1 + bottom



F-90

SPECIFICATIONS

	F-90
Dimensions, H x W x D, mm	652 x 595 x 642
Gross volume, Litres	90
Net volume, Litres	76
Temperature range, C°	-5 to -25
Number of shelves	1 + bottom



ULTF-37

SPECIFICATIONS

	ULTF-37
Dimensions, H x W x D, mm	650 x 630 x 675
Gross volume, Litres	37
Net volume, Litres	37
Temperature range, C°	-65 to -86
Capacity, Cryoboxes "2	33

+20°C



-86°C



LOW NOISE



ENVIRONMENT FRIENDLY



USER FRIENDLY

ISO 9001,
13485, 14001 &
18001 OHSAS
Certified

PROFESSIONAL SOLUTIONS

Vestfrost Solutions is a global developer and manufacturer of innovative and efficient refrigerators and freezers for the professional market. With more than 50 years' experience, we have strong records within manufacturing supplied to the global professional industry, where temperatures must be reliable and precise at all times. We employ numerous initiatives to reduce costs of product usage and to minimize environmental impact.

The continuously expanding product line is based on extensive research and thoroughly tested technology; all products honoring the highest requirements in terms of safety, quality and usability.

For more information, please visit our website:
www.vestfrostsolutions.com



ISO 13485
Management System Certification
BUREAU VERITAS
Certification Denmark A/S

