

PRODUCT

MF 314 VACCINE/ICEPACK FREEZER

The MF series provides vaccine storage and icepack freezing with stable temperature control in the range -18°C to -28°C. The units' thick insulation ensures exceptional hold-over time in case of power cut.

SPECIFICATIONS

Gross volume, litres (cu. ft.)	281 (9.92)
Icepack freezing capacity, kg/24h	7.2
Icepack storage capacity, pcs	256 x 0.6 L
Temperature range (+43°C AMB)	-18° to -28°C
Hold-over time, hours	4.0
Energy consumption - stable, kWh/24h	4.23
Energy consumption - cool down, kWh/24h	4.24
Energy consumption - freezing, kWh/24h	4.24
Refrigerant	R600a
WHO PQS code	E003/023
Climate class	T



FEATURES

Storage baskets, no	3
Temperature control	Automatic
Lock + keys	YES

DIMENSIONS

Height, mm (inches)	850 (33.46)
Width, mm (inches)	1.560 (61.41)
Depth, mm (inches)	600 (23.62)
Gross weight, kg (lbs)	108 (238)
Net weight, kg (lbs)	82 (180)
Shipping volume, m ³	1.16

LOADING QUANTITIES

Qty. per 20' / 40' container	18 / 42
Qty. per truck (EU standard)	48

Products are subject to change due to enhancements and continuous development. Vestfrost Solutions reserve the right to alter any information, without further notice.



COLD CHAIN SOLUTIONS

VESTFROST
SOLUTIONS

A/S Vestfrost
Falkevej 12
DK-6705 Esbjerg Ø
www.vestfrostsolutions.com