



## E003:

## Refrigerators and freezers

<b>PQS code:</b>	<b>E003/070</b>
<b>Type of appliance:</b>	Combined icelined refrigerator/waterpacks freezer
<b>Manufacturer's reference:</b>	VLS 064 RF AC
<b>Manufactured in:</b>	Denmark
<b>Company:</b>	Vestfrost Solutions
<b>Address:</b>	A/S Vestfrost Falkevej 12 DK-6705 Esbjerg Ø Denmark

**Telephone:** +45 7914 2555; +45 7914 2248

**Email:** [clc@vestfrostsolutions.com](mailto:clc@vestfrostsolutions.com);  
[bjn@vestfrostsolutions.com](mailto:bjn@vestfrostsolutions.com)

**Web address:** <http://vestfrostsolutions.com/>

## Specifications

<b>Climate zone:</b>	Hot	<b>Min rated ambient temp:</b>	+10°C
<b>Refrigerant:</b>	R600a	<b>Energy source:</b>	Electric Mains 220-240V 50/60Hz or Electric Mains 115/60Hz optional
<b>Appliance tested at:</b>	+43°C	<b>Ext dimensions (HxLxD)</b>	84.5 x 72.5x 71.0 cm
<b>Performance at:</b>	<b>+43°C</b>	<b>Fuel and cycle type:</b>	Electric - compression

### REFRIGERATOR

### FREEZER

**Freeze protection: Grade A**

**Vaccine storage capacity:** 52.5L  
(Litres)

**Gross volume (Litres):** 75L

**Holdover time (hours):** 20.2

**Energy consumption, stable running (kWh/24 hours):** 1.50

**Energy consumption, cool down test (kWh/24 hours):** 1.51

**Gross volume:** 5.1L

**Waterpack freezing capacity:** 1.6Kg  
6 x 0.6L

**Waterpack storage capacity:**

**Energy consumption during freezing:** 1.80

---

Equipped with unique QR code that can be used with Vestfrost Solutions Unit tracker APP, containing various installation/service features.  
Electrical components positioned in corrosion free metal junction box, protecting for water, dust intrusion, and eases service/maintenance.  
For easy customization to local standard plug, supplied with appliance power in-let + appliance power plug/cable.  
Supplied with adjustable legs, elevating appliance to protect against water/rubbish etc.

---

**Comments:**

**Accessories:**

1 x Vaccine storage top basket; 1 x Bottom grating; 1 x Voltage stabilizer; 1 x 30 day electronic temperature logger; 1 x Holder for temperature logger; Self-regulating temperature control; Electronic temperature reading; 2 x Keys; Instruction manual in 6 UN languages; Technical manual in 6 UN languages; Plywood packing/crate

---

**Spare parts (ref, price):**

Included in the UNICEF Supply Division spare parts kit:  
1 No Compressor - NL7,3MF, Complete - Part No.8-036038987  
1 No Thermostat, XR50CX - unidad de refrigeración - Part No.7095229  
1 No Thermostat, XR50CX - freezer - Part No.7095209  
1 No Starting device, 117U6016 - Part No.6520360  
1 No Thermometer, solar - Part No.7020287-03

Other parts are available as spares on request.

---

**Shipping volume:**

0,58 m3

**Shipping weight:**

98 kg

**Price year**

2016

**Incoterms**

FCA

**Price**

---

**Quality standard:**

- ISO 9001:2008 - ISO 14001

**Specification ref.**

-

---

**Verification report:**

300-KLAB-15-181

**Verification laboratory:**

Danish Technological Institute

---

**Current PQS status:**

*pre-qualified: 08 Apr 2016*

**Validity until:** *May, 2021*

---

Note: If Current PQS status is 'Suspended' or 'Withdrawn', this product is NOT to be purchased.