

# Datasheet



## BIOMEDICAL

### CATEGORY - FREEZER

### MODEL - VT 306

#### Controller

Defrost	Manual
High/Low Temp. Alarm	Yes
Open Door Alarm	No
Probe Failure Alarm	Yes
Power Failure Alarm	No
USB Connection	No
Temperature Graph	No

#### Features

Lock	Yes
LED Light	No
Battery Backup for Controller, 48h	No
Porthole	Yes
Dry Contact	No
Compressor	1

#### Construction

Outer Dimensions, HxWxD	885x1265x605
Inner Dimensions, HxWxD	650x1100x437
Volume Gross/net, L	296 / 284
Net Weight, kg	63,4
Material Inner Cabinet	Painted Steel
Material Outer Cabinet	Painted Steel
Insulation Thickness, mm	80
Insulation Type	Polyurethane with Cyclopentane
Refrigerant	R290
Refrigeration System	Static
Feet / Castors	Castors



LOW NOISE



ENVIRONMENT FRIENDLY



USER FRIENDLY

#### Operation

Temperature Range, C°	-25 to -45
Noise, dB	55
Energy Consumption 24/h, kWh	3,7
Frequency, Hz	50
Power Supply, V	220 / 240

#### Storage

Number of Baskets (Standard/max)	1 / 1
Shelf material	Steel coated with Plastic Powder
Max Loading pr. Basket, kg	16
Max no. Drawers Available	N/A
Drawer Material	N/A
Max Loading pr. Drawer, kg	N/A