

BETWEEN
-86°C

-60°C

BIOMEDICAL SOLUTIONS
**NEW ULTRA LOW
TEMPERATURE FREEZER**

NEW ULT FREEZER

ULTF-37

Introducing our new ULT Freezer; ULTF-37. Built as a Tabletop it can store your samples right next to your workstation.

With a small footprint the "personal freezer" fits perfectly on the counter where you work. This reduces the time you have to spend on retrieving and storing samples from a distant, centralized location. It also contributes positively to the integrity of the samples, by reducing the time they are exposed to room temperature, as you can easily and quickly store them again once you are done testing.

Assuming performance and features are similar, much of the perceived value in a ULT purchase decision, resides around the ratio between acquisition price and volume. However, if you want just the right ULT freezer for your needs, the picture is much more nuanced.

Many laboratories have limited space and with a multitude of necessary equipment, more is, not always better, when evaluating your freezer options. It is also relevant to ask yourself how much capacity you really need for your daily work with samples. This is due to the fact that a freezer will perform better in terms of temperature conditions, and in energy consumption, when it is almost at full capacity.



THE SPECIFICATIONS

ULTF-37

37 LITRES GROSS VOLUME

- Vacuum Insulated Panels
- Variable Speed Compressor
- Hot wire in front panel

H: 650 mm W: 630 mm D: 675 mm



LOW NOISE
<50 dB



ENVIRONMENT FRIENDLY
Natural Refrigerant
Approx. 3,2 kWh/24h



USER FRIENDLY
Plug & Play
Ergonomic door handle
Filterless construction
Adjustable legs for leveling

TECHNICAL SPECIFICATIONS	ULTF-37
Dimensions, H x W x D, mm	650 x 630 x 675
Gross volume, Litres	37
Net volume, Litres	37
Temperature range, C°	-60 to -86

DESIGNED FOR FREQUENT USE

The size and dimensions of the ULTF-37 is ideal for reducing a rise in temperature associated with frequent door openings. Furthermore we have installed a hot wire in our panels to combat the icing that inevitably propagate as a result of the door being opened.



Personal Freezer

Store & retrieve the samples you work with daily, easily and quickly, right at your desk

Impressive capacity

The tabletop has an impressive capacity of 37L

Net volume translating to:

- 33x 2" Cryoboxes
- Up to 3.300 2ml vials

RELIABILITY IS STILL THE KEY

None of the beforementioned features matter if the product doesn't deliver on its most basic premise - Safe and reliable storage. We have not forgotten this purpose, which is why all of our units undergo a 24h test before they are shipped.

And when you plug in your ULTF-37, it will feature all the necessary alarms to ensure that you are alerted if temperatures are less than ideal in your freezer.

Controller Features

- Power Failure alarm
- Probe Failure alarm
- High/low Temperature alarm
- Door Opening alarm
- Data Logging with USB Extraction



MADE IN
DENMARK