

VT 306

LOW TEMPERATURE CHEST FREEZER

The VT low temperature freezers creates the possibility to maintain temperatures as low as -60°C. Supreme stability, reliability, user-friendliness and ease of cleaning make these freezers an ideal solution.



-10°C
-45°C

Products are subject to change due to enhancements and continuous development. Vestfrost Solutions reserve the right to alter any information, without further notice.

| CONSTRUCTION | | |
|---|--|---------------------|
| Outer Dimensions, HxWxD, mm | | 885x1265x605 |
| Inner Dimensions, HxWxD, mm | | 635x1100x440 |
| PACKAGING | | |
| Dimensions HxWxD, mm | | - |
| VOLUME | | |
| Volume, Gross, L | | 296 |
| Volume, Net, L | | 296 |
| WEIGHT | | |
| Weight, Gross, kg | | 77 |
| Weight, Net, kg | | 71 |
| BASKETS | | |
| No. of baskets, pcs | | 1 |
| TEMPERATURE | | |
| Ambient temperature, °C | | +10 to +30 |
| Climate class | | - |
| Temperature range, °C | | -10 to -45 |
| COOLING SYSTEM | | |
| Refrigerant, type | | R290 |
| Refrigerant amount, gr | | - |
| Static/ventilated | | Static |
| ELECTRICAL | | |
| Voltage, V | | 220-240 |
| Frequency, Hz | | 50 |
| Voltage 110/60, V-Hz | | Option |
| Type of thermometer | | Electronic |
| Type of thermostat | | Electronic |
| Potential free contact for remote control | | Yes |
| Porthole for sensor, ø12 mm | | Yes |
| ENERGY | | |
| Energy consumption, kWh/24h | | 3,74 |
| Energy consumption, kWh/year | | - |
| Energy class | | not needed |
| EEl | | not needed |
| Input power, watt | | 415 |
| FEATURES | | |
| Cabinet color | | White |
| Inside material/color | | White painted steel |
| Door or lid color/type | | White, lift up |
| Alternative lid | | n/a |
| Lock | | Yes |
| Castors/feet | | Castors |
| Defrost & removal of water | | Manual |
| Insulation, mm | | 80 |
| OPTION | | |
| Graphics | | Option |
| Baskets | | - |
| SHIPPING | | |
| Quantity in standard truck | | 120 |
| Quantity in 20 ft container | | - |
| Quantity in 40 ft HC container | | 54 |