

# VT 307

## LOW TEMPERATURE CHEST FREEZER

The VT low temperature freezers creates the possibility to maintain temperatures as low as -60°C. Supreme stability, reliability, user-friendliness and ease of cleaning make these freezers an ideal solution.



Products are subject to change due to enhancements and continuous development. Vestfrost Solutions reserve the right to alter any information, without further notice.

| CONSTRUCTION                              |                     |
|-------------------------------------------|---------------------|
| Outer Dimensions, WxDxH, mm               | 1260x600x885        |
| Inner Dimensions, WxDxH, mm               | 1100x440x635        |
| PACKAGING                                 |                     |
| Dimensions WxDxH, mm                      | -                   |
| VOLUME                                    |                     |
| Volume, Gross/Net, L                      | 296/284             |
| WEIGHT                                    |                     |
| Weight, Gross/Net, kg                     | 76/71               |
| BASKETS                                   |                     |
| No. of baskets, pcs                       | 3                   |
| TEMPERATURE                               |                     |
| Ambient temperature, °C                   | +10 to +30          |
| Climate class                             | -                   |
| Temperature range, °C                     | -40 to -60          |
| COOLING SYSTEM                            |                     |
| Refrigerant, type                         | Mix of hydrocarbons |
| Refrigerant amount, gr                    | -                   |
| Static/ventilated                         | Static              |
| ELECTRICAL                                |                     |
| Voltage, V                                | 220-240             |
| Frequency, Hz                             | 50                  |
| Voltage 110/60, V-Hz                      | Option              |
| Variable speed compressor                 | option              |
| Type of controller                        | Dixell XR30CX       |
| Potential free contact for remote control | Yes                 |
| Porthole for sensor, ø12 mm               | Yes                 |
| ENERGY                                    |                     |
| Energy consumption, kWh/24h               | 5,56                |
| Energy cons. 24 h variable speed          | 3,68                |
| Input power, watt                         | -                   |
| FEATURES                                  |                     |
| Cabinet color                             | White               |
| Inside material/color                     | White painted steel |
| Door or lid color/type                    | White, lift up      |
| Alternative lid                           | n/a                 |
| Lock                                      | Yes                 |
| Castors/feet                              | Castors             |
| Defrost & removal of water                | Manual              |
| Insulation, mm                            | 80                  |
| OPTION                                    |                     |
| Controller language                       | 3 digits            |
| USB connection                            | No                  |
| Logging                                   | No                  |
| Temperature graph                         | No                  |
| High / low temp. alarm                    | Yes                 |
| Probe failure alarm                       | Yes                 |
| Power failure alarm                       | No                  |
| SHIPPING                                  |                     |
| Quantity in standard truck                | 120                 |
| Quantity in 20 ft container               | -                   |
| Quantity in 40 ft HC container            | 54                  |