

# PRODUCT

## MF 314 VACCINE/ICEPACK FREEZER

The MF series provides vaccine storage and icepack freezing with stable temperature control in the range -15°C to -25°C. The units' thick insulation ensures exceptional hold-over time in case of power cut.



### SPECIFICATIONS

Gross volume	281
Icepack freezing capacity, kg/24h	7.2
Icepack storage capacity, pcs	256 x 0.6L
Temperature range (+43C AMB)	-15 to -25C
Hold-over time, hours	4.0
Energy consumption - stable, kWh/24h	4.23
Energy consumption - cool down, kWh/24h	4.29
Refrigerant	R600a
WHO PQS code	E003/023
Climate class	T

### FEATURES

Storage baskets, no	3
Temperature control	Automatic
Lock + keys	YES
Adjustable legs	YES, 10cm
Junction box	YES
Safety thermostat	YES
Appliance plugs	M,G,F
Industrial grade locking mechanism	YES
Fridge-tag 2E	YES

### DIMENSIONS

Height, mm	910
Width, mm	1560
Depth, mm with handle	700
Gross weight, kg	106
Net weight, kg	80
Shipping volume m3	1.21

### LOADING QUANTITIES

Qty. per 20' / 40' container	18 / 42
Qty. per truck (EU standard)	48

### VOLTAGE REGULATION RANGE

220-240V/50/60Hz	110-278V/50/60Hz
100-127V/50/60Hz	82-159V/50/60Hz



Products are subject to change due to enhancements and continuous development. Vestfrost Solutions reserve the right to alter any information, without further notice.



**COLD CHAIN SOLUTIONS**

**VESTFROST**  
SOLUTIONS

A/S Vestfrost  
Falkevej 12  
DK-6705 Esbjerg Ø  
www.vestfrostsolutions.com