

Refrigerants

Overview of Refrigerant quantity in each model, as well as their respective Global Warming Potential (GWP) & Ozone Depletion Potential (ODP).



ENVIRONMENT FRIENDLY

REFRIGERATORS	Type	Refrigerant	Qty (Gram)	GWP	ODP
R-90	Tabletop	R600a	39	3	0
R-138	Undercounter	R600a	47	3	0
157	Undercounter	R600a	40	3	0
337	Vertical	R600a	50	3	0
397	Vertical	R600a	60	3	0
427	Vertical	R600a	60	3	0
FREEZERS	Type	Refrigerant	Qty (Gram)	GWP	ODP
F-90	Tabletop	R600a	40	3	0
252	Vertical	R290	70	3,3	0
146	Chest	R290	50	3,3	0
206	Chest	R290	60	3,3	0
306	Chest	R290	65	3,3	0
406	Chest	R290	70	3,3	0
546	Chest	R290	120	3,3	0
254	Vertical	Nature R	149	4,4	0
256	Vertical	Nature R	167	4,4	0
75	Chest	Nature R	100	4,4	0
147	Chest	Nature R	137	4,4	0
207	Chest	Nature R	88	4,4	0
307	Chest	Nature R	166	4,4	0
407	Chest	Nature R	195	4,4	0
547	Chest	Nature R	170	4,4	0
ULT	Type	Refrigerant	Qty (Gram)	GWP	ODP
ULTF-37	Vertical	Nature R	96	4,4	0
098	Vertical	Nature R	109	4,4	0
258	Vertical	Nature R	172	4,4	0
78	Vertical	Nature R	121	4,4	0
208	Vertical	Nature R	90	4,4	0
308	Vertical	Nature R	143	4,4	0