

## Refrigerants

Overview of Refrigerant quantity in each model, as well as their respective Global Warming Potential (GWP) & Ozone Depletion Potential (ODP).



ENVIRONMENT FRIENDLY

### Refrigerators +2°C to +20°C

MODEL	TYPE	REFRIGERANT	QTY (GRAM)	GWP	ODP
R-90i	Tabletop	R600a	39	3	0
R138i	Undercounter	R600a	47	3	0
R-427i	Upright	R600a	60	3	0

### Freezers -5°C to -25°C

MODEL	TYPE	REFRIGERANT	QTY (GRAM)	GWP	ODP
F-90i	Tabletop	R600a	38	3	0
F138i	Undercounter	R600a	43	3	0
F-427i	Upright	R290	80	3	0

### ULT Freezers -20°C to -86°C

MODEL	TYPE	REFRIGERANT	QTY (GRAM)	GWP	ODP
ULTF-37i	Tabletop	Nature R 2	96	4,4	0
VTS 098	Undercounter	Nature R 2	109	4,4	0
ULTF-C74i	Chest	Nature R 2	121	4,4	0
ULTF-C198i	Chest	Nature R 2	90	4,4	0
ULTF-C296i	Chest	Nature R 2	143	4,4	0
ULTF-C383i	Chest	Nature R 2	161	4,4	0