

VT 75

LOW TEMPERATURE CHEST FREEZER

The VT low temperature freezers creates the possibility to maintain temperatures as low as -60°C. Supreme stability, reliability, user-friendliness and ease of cleaning make these freezers an ideal solution.



Products are subject to change due to enhancements and continuous development. Vestfrost Solutions reserve the right to alter any information, without further notice.

CONSTRUCTION		
Outer Dimensions, WxDxH, inch		21.85×21.85×34.25
Inner Dimensions, WxDxH, inch		15.35×15.35×18.31
PACKAGING		
Dimensions WxDxH, inch		-
VOLUME		
Volume, Gross, cu ft		2.51 ft³
Volume, Net, cu ft		2.51 ft³
WEIGHT		
Weight, Gross, lb		105.82 lb
Weight, Net, lb		92.59 lb
BASKETS		
No. of baskets, pcs		1
TEMPERATURE		
Ambient temperature, °F		50 °F to 86 °F
Climate class		-
Temperature range, °F		40 °F to -76 °F
COOLING SYSTEM		
Refrigerant, type		Mix of hydrocarbons
Refrigerant amount, gr		-
Static/ventilated		Static
ELECTRICAL		
Voltage, V		110 Hz
Frequency, Hz		60 V
Voltage 110/60, V-Hz		Option
Type of thermometer		Electronic
Type of thermostat		Electronic
Potential free contact for remote control		Yes
Porthole for sensor, ø12 mm		Yes
ENERGY		
Energy consumption, kWh/24h		2,2
Energy consumption, kWh/year		-
Energy class		not needed
EEL		not needed
Input power, watt		415
FEATURES		
Cabinet color		White
Inside material/color		White painted steel
Door or lid color/type		White, lift up
Alternative lid		n/a
Lock		Yes
Castors/feet		Feet
Defrost & removal of water		Manual
Insulation, mm		80
OPTION		
Graphics		Option
Baskets - max. qty		-
SHIPPING		
Quantity in standard truck		288
Quantity in 20 ft container		-
Quantity in 40 ft HC container		114