

# VT 406

## LOW TEMPERATURE CHEST FREEZER

The VT low temperature freezers creates the possibility to maintain temperatures as low as -60°C. Supreme stability, reliability, user-friendliness and ease of cleaning make these freezers an ideal solution.



-10°C  
-45°C

Products are subject to change due to enhancements and continuous development. Vestfrost Solutions reserve the right to alter any information, without further notice.

| CONSTRUCTION                              |  |                     |
|-------------------------------------------|--|---------------------|
| Outer Dimensions, WxDxH, mm               |  | 1562x605x885        |
| Inner Dimensions, WxDxH, mm               |  | 1400x440x635        |
| PACKAGING                                 |  |                     |
| Dimensions WxDxH, mm                      |  | -                   |
| VOLUME                                    |  |                     |
| Volume, Gross, L                          |  | 384                 |
| Volume, Net, L                            |  | 384                 |
| WEIGHT                                    |  |                     |
| Weight, Gross, kg                         |  | 83                  |
| Weight, Net, kg                           |  | 77                  |
| BASKETS                                   |  |                     |
| No. of baskets, pcs                       |  | 1                   |
| TEMPERATURE                               |  |                     |
| Ambient temperature, °C                   |  | +10 to +30          |
| Climate class                             |  | -                   |
| Temperature range, °C                     |  | -10 to -45          |
| COOLING SYSTEM                            |  |                     |
| Refrigerant, type                         |  | R290                |
| Refrigerant amount, gr                    |  | -                   |
| Static/ventilated                         |  | Static              |
| ELECTRICAL                                |  |                     |
| Voltage, V                                |  | 220-240             |
| Frequency, Hz                             |  | 50                  |
| Voltage 110/60, V-Hz                      |  | Option              |
| Type of thermometer                       |  | Electronic          |
| Type of thermostat                        |  | Electronic          |
| Potential free contact for remote control |  | Yes                 |
| Porthole for sensor, ø12 mm               |  | Yes                 |
| ENERGY                                    |  |                     |
| Energy consumption, kWh/24h               |  | 4,09                |
| Energy consumption, kWh/year              |  | -                   |
| Energy class                              |  | not needed          |
| EEl                                       |  | not needed          |
| Input power, watt                         |  | 415                 |
| FEATURES                                  |  |                     |
| Cabinet color                             |  | White               |
| Inside material/color                     |  | White painted steel |
| Door or lid color/type                    |  | White, lift up      |
| Alternative lid                           |  | n/a                 |
| Lock                                      |  | Yes                 |
| Castors/feet                              |  | Castors             |
| Defrost & removal of water                |  | Manual              |
| Insulation, mm                            |  | 80                  |
| OPTION                                    |  |                     |
| Graphics                                  |  | Option              |
| Baskets - max. qty                        |  | -                   |
| SHIPPING                                  |  |                     |
| Quantity in standard truck                |  | 96                  |
| Quantity in 20 ft container               |  | -                   |
| Quantity in 40 ft HC container            |  | 42                  |