



E003:

Refrigerators and freezers

PQS code:	E003/024
Type of appliance:	Vaccine/waterpacks freezer
Manufacturer's reference:	MF 114
Manufactured in:	Denmark
Company:	Vestfrost Solutions
Address:	A/S Vestfrost Falkevej 12 DK-6705 Esbjerg Ø Denmark
Telephone:	+45 7914 2555; +45 7914 2248
Email:	clc@vestfrostsolutions.com ; bjn@vestfrostsolutions.com
Web address:	http://vestfrostsolutions.com/

Specifications

Climate zone:	Hot	Min rated ambient temp:	N/A
Refrigerant:	R600a	Energy source:	Electric Mains 220-240V 50/60Hz or Electric Mains 115/60Hz optional
Appliance tested at:	+43°C	Ext dimensions (HxLxD)	84 x 72 x 70 cm
Performance at:	+43°C	Fuel and cycle type:	Electric - compression

REFRIGERATOR

FREEZER

Freeze protection: N/A

Vaccine storage capacity: (Litres)	N/A		
Gross volume (Litres):	N/A	Gross volume:	105 L
Holdover time (hours):	2.8 hrs	Waterpack freezing capacity:	7.2 kg/24h 64 x 0.6 L
Energy consumption, stable running (kWh/24 hours):	2.24 kWh/24hrs	Waterpack storage capacity:	
Energy consumption, cool down test (kWh/24 hours):	2.41 kWh/24hrs	Energy consumption during freezing:	3.33 kWh/24hrs

