

# VT 206

## LOW TEMPERATURE CHEST FREEZER

The VT low temperature freezers creates the possibility to maintain temperatures as low as -60°C. Supreme stability, reliability, user-friendliness and ease of cleaning make these freezers an ideal solution.



Products are subject to change due to enhancements and continuous development. Vestfrost Solutions reserve the right to alter any information without further notice.

| CONSTRUCTION                              |                     |
|---|---------------------|
| Outer Dimensions, HxWxD, mm               | 885x923x605         |
| Inner Dimensions, HxWxD, mm               | 635x760x440         |
| PACKAGING                                 |                     |
| Dimensions HxWxD, mm                      | -                   |
| VOLUME                                    |                     |
| Volume, Gross, L                          | 198                 |
| Volume, Net, L                            | 198                 |
| WEIGHT                                    |                     |
| Weight, Gross, kg                         | 63                  |
| Weight, Net, kg                           | 58                  |
| BASKETS                                   |                     |
| No. of baskets, pcs                       | 1                   |
| TEMPERATURE                               |                     |
| Ambient temperature, °C                   | +10 to +35          |
| Climate class                             | -                   |
| Temperature range, °C                     | -10 to -45          |
| COOLING SYSTEM                            |                     |
| Refrigerant, type                         | R290                |
| Refrigerant amount, gr                    | -                   |
| Static/ventilated                         | Static              |
| ELECTRICAL                                |                     |
| Voltage, V                                | 220-240             |
| Frequency, Hz                             | 50                  |
| Voltage 110/60, V-Hz                      | Option              |
| Type of thermometer                       | Electronic          |
| Type of thermostat                        | Electronic          |
| Potential free contact for remote control | Yes                 |
| Porthole for sensor, ø12 mm               | Yes                 |
| ENERGY                                    |                     |
| Energy consumption, kWh/24h               | 3,14                |
| Energy consumption, kWh/year              | -                   |
| Energy class                              | not needed          |
| EEL                                       | not needed          |
| Input power, watt                         | 415                 |
| FEATURES                                  |                     |
| Cabinet color                             | White               |
| Inside material/color                     | White painted steel |
| Door or lid color/type                    | White, lift up      |
| Alternative lid                           | n/a                 |
| Lock                                      | Yes                 |
| Castors/feet                              | Castors             |
| Defrost & removal of water                | Manual              |
| Insulation, mm                            | 80                  |
| OPTION                                    |                     |
| Graphics                                  | Option              |
| Baskets                                   | -                   |
| SHIPPING                                  |                     |
| Quantity in standard truck                | 168                 |
| Quantity in 20 ft container               | -                   |
| Quantity in 40 ft HC container            | 78                  |