

VT 307

LOW TEMPERATURE CHEST FREEZER

The VT low temperature freezers creates the possibility to maintain temperatures as low as -60°C. Supreme stability, reliability, user-friendliness and ease of cleaning make these freezers an ideal solution.



+40°F
-76°F

Products are subject to change due to enhancements and continuous development. Vestfrost Solutions reserve the right to alter any information, without further notice.

CONSTRUCTION		
Outer Dimensions, WxDxH, inch		49.80×23.82×34.84
Inner Dimensions, WxDxH, inch		43.31×17.32×25.00
PACKAGING		
Dimensions WxDxH, inch		-
VOLUME		
Volume, Gross, cu ft		10.46 ft³
Volume, Net, cu ft		10.46 ft³
WEIGHT		
Weight, Gross, lb		169.75 lb
Weight, Net, lb		156.53 lb
BASKETS		
No. of baskets, pcs		3
TEMPERATURE		
Ambient temperature, °F		50 °F to 86 °F
Climate class		-
Temperature range, °F		40 °F to -76 °F
COOLING SYSTEM		
Refrigerant, type		Mix of hydrocarbons
Refrigerant amount, gr		-
Static/ventilated		Static
ELECTRICAL		
Voltage, V		110 Hz
Frequency, Hz		60 V
Voltage 110/60, V-Hz		Option
Type of thermometer		Electronic
Type of thermostat		Electronic
Potential free contact for remote control		Yes
Porthole for sensor, ø12 mm		Yes
ENERGY		
Energy consumption, kWh/24h		5,56
Energy consumption, kWh/year		-
Energy class		not needed
EEL		not needed
Input power, watt		-
FEATURES		
Cabinet color		White
Inside material/color		White painted steel
Door or lid color/type		White, lift up
Alternative lid		n/a
Lock		Yes
Castors/feet		Castors
Defrost & removal of water		Manual
Insulation, mm		80
OPTION		
Graphics		Option
Baskets - max. qty		-
SHIPPING		
Quantity in standard truck		120
Quantity in 20 ft container		-
Quantity in 40 ft HC container		54