

VT 407

LOW TEMPERATURE CHEST FREEZER

The VT low temperature freezers creates the possibility to maintain temperatures as low as -60°C. Supreme stability, reliability, user-friendliness and ease of cleaning make these freezers an ideal solution.



+40°F
-76°F

Products are subject to change due to enhancements and continuous development. Vestfrost Solutions reserve the right to alter any information, without further notice.

CONSTRUCTION	
Outer Dimensions, WxDxH, inch	61.50×23.82×34.84
Inner Dimensions, WxDxH, inch	55.12×17.32×25.00
PACKAGING	
Dimensions WxDxH, inch	-
VOLUME	
Volume, Gross, cu ft	13.55 ft³
Volume, Net, cu ft	13.55 ft³
WEIGHT	
Weight, Gross, lb	182.98 lb
Weight, Net, lb	169.75 lb
BASKETS	
No. of baskets, pcs	4
TEMPERATURE	
Ambient temperature, °F	50 °F to 86 °F
Climate class	-
Temperature range, °F	40 °F to -76 °F
COOLING SYSTEM	
Refrigerant, type	Mix of hydrocarbons
Refrigerant amount, gr	-
Static/ventilated	Static
ELECTRICAL	
Voltage, V	110 Hz
Frequency, Hz	60 V
Voltage 110/60, V-Hz	Option
Type of thermometer	Electronic
Type of thermostat	Electronic
Potential free contact for remote control	Yes
Porthole for sensor, ø12 mm	Yes
ENERGY	
Energy consumption, kWh/24h	7,08
Energy consumption, kWh/year	-
Energy class	not needed
EEL	not needed
Input power, watt	-
FEATURES	
Cabinet color	White
Inside material/color	White painted steel
Door or lid color/type	White, lift up
Alternative lid	n/a
Lock	Yes
Castors/feet	Castors
Defrost & removal of water	Manual
Insulation, mm	80
OPTION	
Graphics	Option
Baskets - max. qty	-
SHIPPING	
Quantity in standard truck	96
Quantity in 20 ft container	-
Quantity in 40 ft HC container	42